

# wish Bone E Mg

# MAGNESIUM-BASED BONE VOID FILLER

# **MIXING INSTRUCTIONS**



# **TABLE OF CONTENTS**

## Introduction

Product Overview	З
Kit Components	З
Features & Benefits	4
Why Magnesium	5
Indications	5
Product Comparison	6

# **Mixing Instructions**

Moldable	7
Injectable	-9
Working Guidelines - Curing Time	10

# **Order Information**

WishBONE Mg BVF Kits (3)	11
--------------------------	----

# MAGNESIUM-BASED BONE VOID EILLER

# Injectable and Moldable BVF | Full Kits for Mixing + Application

### High Strength. Lightweight. Resorbable.

WishBONE Mg is a moldable/injectable Magnesium-based bone void filler that has a unique resorption profile that provides stability while also increasing cell proliferation, advancement of mineralization with a result of enhanced bone regeneration for multiple types of orthopedic applications.<sup>1-8, 13-15</sup>

WishBONE Mg made from a pre-measured blend of magnesium, phosphates and a pre-measured proprietary solution. When mixed and molded/injected according to the instructions for use, the product will harden in situ at the defect site.<sup>12</sup>

- 80% resorbable in 26 weeks<sup>9</sup>
- Remodels to normal bone<sup>10, 11</sup>
- Quicker time to union compared to calcium-based BVFs<sup>9, 12</sup>



# FEATURES AND BENEFITS

### **INJECTABLE & MOLDABLE**

Designed with optimized consistency for use in multiple orthopaedic applications<sup>13, 15</sup>

### ENHANCED BONE REGENERATION

Greater than 80% bone remodeling in 26 weeks9

### FULLY SYNTHETIC MATERIAL

Enhanced quality control and increased product availability with reduced product morbidity compared to human tissue

### RADIOPAQUE

Product easily identifiable in situ

### OSTEO-CONDUCTIVE/HIGH COMPRESSIVE STRENGTH

Surface topography to support bone formation enhancing structural stability and biocompatibility<sup>12, 14</sup>

### **TEMPERATURE SETTING CONTROL**

Designed for optimal workability and curing time<sup>13, 14</sup>

### **THIXOTROPIC PROPERTIES**

Easily manageable curing time designed to allow for intra-operative flexibility and reduction of waste<sup>13</sup>

### **EXCELLENT BINDING CHARACTERISTICS**

Optimal product stability and fixation at operative site<sup>9, 12</sup>

**CASE STUDY** – 2 weeks, 6 weeks, 3 months and 6 months Tibia Plateau Fracture, treated with 10cc OsteoCrete



2 weeks

6 weeks



3 months

6 months

Images courtesy of Dr. Wetzel – University Hospital – Cleveland, OH

Claims based on critically sized rabbit lateral condyle defect model, rabbit anterior cruciate ligament reconstruction, equine metacarpal and metatarsal fracture fixation, and equine metatarsal osteotomy. It is unknown how results from the rabbit or equine models compare with clinical results in humans.



# WHY MAGNESIUM?

A natural nutrient essential for building healthy bones, magnesium is an attractive orthopedic alternative because, unlike stainless steel or titanium, it is biodegradable – eliminating the need for an invasive procedure to remove surgical hardware after initial surgery.

- Plays a role in the active transport of calcium and potassium ions across cell membranes<sup>1</sup>
- Contributes to the structural development of bone<sup>1-8</sup>
- Approximately 60% of Mg in the body resides in bones<sup>16</sup>



For product information, including indications, contraindications, warnings, precautions and potential adverse effects, visit WishBone Medical's Instructions for Use page online: www.WishBoneMedical.com/IFU.

# **PRODUCT COMPARISON:** WISHBONE MG VS. COMPETITOR OFFERINGS

Product	Synthetic	Injectable	Moldable	Radiopaque	Osteo- Conductive	High Compressive Strength	Mixing Delivery System Provided	Temperature Setting Control	Binding Characteristics	Thixotropic Waste Reduction Properties	Enhanced Bone Regeneration
WishBONE Mg (WishBone Medical)	1	1	1	1	1	1	1	1	1	✓	1
CERAMENT™ (BoneSupport)'⁰	1	1	1		1		1				
Pro-Dense™ (Wright Medical)²⁰	1	1	1	1	1	1	1				
HydroSet™ (Stryker)²¹	1	1	1	1	1	1	1				
Norian <sup>®</sup> (DePuy Synthes) <sup>22</sup>	1	1	1	1	1	1		1			
AccuFiLL® BSM (Zimmer Biomet) <sup>23</sup>	1	1		1	1		1	1			
Pro-Dense™ (Wright Medical)²⁰	1	1	1	1	1	1	1	1			
Callos® Inject (Acumed/Skeletal Kinetics) <sup>24</sup>	1	1	1	1	1	1	1	1			

# MIXING INSTRUCTIONS: MOLDABLE

**Preparation:** The surgical field should be irrigated to remove any loose debris and dried prior to placement of WishBONE Mg. Before mixing, ensure that WishBONE Mg is equilibrated to room temperature: (18-23°C / 65-73°F)

 Open the wishBONE Mg powder pouch and premeasured liquid solution and pour both into the sterile basin (18-23°C / 65-73°F).



With the spatula, mix until powder and liquid solution for 2 minutes until product is thoroughly mixed into a smooth material (18-23°C / 65-73°F).



- Wait until a skin is formed on the product during the curing process (18-23°C / 65-73°F).
  - **Technique Tip:** Increased ambient temperature of the operating room will accelerate curing time.



After the waiting time, the product will be in a moldable putty form and ready for implantation. It can be contoured manually or with an instrument if desired (18-23°C / 65-73°F).

For mixing and curing times, reference chart on pg. 10.





# MIXING INSTRUCTIONS: INJECTABLE

Remove the white cap assembly at the distal end of the syringe by turning counter clockwise.
Pull the plunger all the way back.



Attach the funnel to the distal end by turning clockwise. Transfer the powder into the syringe. Make sure all powder is transferred.

This might require light tapping of the funnel to completely transfer the material into the syringe.

Add the premeasured liquid solution into the syringe.

Remove the funnel and reattach the white cap assembly. Be sure the cap and wing nut are both secured tightly.





# MIXING INSTRUCTIONS: INJECTABLE (CONT.)

Remove the white cap assembly at the distal end of the syringe by turning counter clockwise. Pull the plunger all the way back. Remove the stopper to enable the mixing element.

**Reminder:** Do not discard the stopper at this point. Mix back and forth for 2 minutes with a twisting motion until powder and liquid solution are thoroughly mixed to a smooth material.

**Technique Tip:** Increased ambient temperature of the operating room will accelerate the curing time.

**Technique Tip:** Continue mixing product to prolong set time.





Pull the mixing element all the way back and reattach the stopper (see step 4). Remove the wing nut from the distal cap and plunge the material to remove the air from syringe. Attach a luer compatible device a begin injecting material as needed.





# WISHBONE MG WORKING GUIDELINES<sup>25</sup>





# **ORDER INFORMATION**

KIT NUMBER	DESCRIPTION
PK-BVF05	WishBONE Mg Bone Void Filler Kit, 5cc
PK-BVF10	WishBONE Mg Bone Void Filler Kit, 10cc
PK-BVF15	WishBONE Mg Bone Void Filler Kit, 15cc

Kit Contents:. Liquid & Powder Set, Mixing Syringe, Funnel, Basin, Spatula & 11GA Cannula



**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

## REFERENCES

- 1. Díaz-Tocados JM, Herencia C, Martínez-Moreno JM, et al. Magnesium Chloride promotes Osteogenesis through Notch signaling activation and expansion of Mesenchymal Stem Cells. Sci Rep. 2017;7(1):7839.
- 2. Yoshizawa S, Brown A, Barchowsky A, Sfeir C. Magnesium ion stimulation of bone marrow stromal cells enhances osteogenic activity, simulating the effect of magnesium alloy degradation. Acta Biomater. 2014;10(6):2834-2842.
- 3. He LY, Zhang XM, Liu B, Tian Y, Ma WH. Effect of magnesium ion on human osteoblast activity. Braz J Med Biol Res. 2016;49(7).
- 4. Wong H, Chu P, Leung K, Cheung M, Luk K, Yeung K. Engineered polycaprolactone-magnesium hybrid biodegradable porous scaffold for bone tissue engineering. Progress in Natural Science: Materials International. 2014;24(5):561-567.
- 5. Jia J, Zhou H, Wei J, et al. Development of magnesium calcium phosphate biocement for bone regeneration. J R Soc Interface. 2010;7(49):1171-1180.
- Zeng D, Xia L, Zhang W, et al. Maxillary sinus floor elevation using a tissue-engineered bone with calciummagnesium phosphate cement and bone marrow stromal cells in rabbits. Tissue Eng Part A. 2012;18(7-8):870881.
- Wu F, Wei J, Guo H, Chen F, Hong H, Liu C. Self-setting bioactive calcium-magnesium phosphate cement with high strength and degradability for bone regeneration. Acta Biomater. 2008;4(6):1873-1884.
- 8. Zhang Z, Yang Z, Chen Z, et al. A study on bone cement containing magnesium potassium phosphate for bone repair. Cogent Biology. 2018;4(1):1-11.
- 9. Ref. Bone Solutions Test Report: 10188 Rev. A / Data on file
- 10. Bertone A, DeMaria M, Johnson A, Weisbrode S, Kowaleski M. Degradable magnesium based cement adheres stainless steel screws into bone. Orthopaedic Research Society; 2006; Chicago, IL.
- 11. Bertone A, Hackett B, Litsky A, Johnson A, Kaeding C, Lally T. A magnesium injectable formulation adheres bone to bone and tendon to bone. Orthopaedic Research Society; 2005; Washington, D.C.
- 12. Waselau M, Samii VF, Weisbrode SE, Litsky AS, Bertone AL. Effects of a magnesium adhesive cement on bone stability and healing following a metatarsal osteotomy in horses. Am J Vet Res. 2007;68(4):370-378.
- 13. Hirvinen LJ, Litsky AS, Samii VF, Weisbrode SE, Bertone AL. Influence of bone cements on bone-screw interfaces in the third metacarpal and third metatarsal bones of horses. Am J Vet Res. 2009;70(8):964-972.
- 14. Kim MS, Kovacevic D, Milks RA, et al. Bone Graft Substitute Provides Metaphyseal Fixation for a Stemless Humeral Implant. Orthopedics. 2015;38(7):e597-603. 13-15
- 15. Gulotta LV, Kovacevic D, Ying L, Ehteshami JR, Montgomery S, Rodeo SA. Augmentation of tendon-to-bone healing with a magnesium-based bone adhesive. Am J Sports Med. 2008;36(7):1290-1297.
- 16. Gröber U, Schmidt J, Kisters K. Magnesium in Prevention and Therapy. Nutrients. 2015;7(9):8199–8226. Published 2015 Sep 23. doi:10.3390/nu7095388
- 17. Yoshizawa et al. Magnesium ion stimulation of bone marrow stromal cells enhances osteogenic activity, stimulating the effect of magneisum alloy degradation. Acta Biomater. 2014; 10(6): 2834-42.
- 18. Wong et al. Engineered polycaprolactone-magnesium hybrid biodegradable porous scaffold for bone tissue engineering. Materials International. 2014; 24: 561-567.
- 19. Ref. Biomet Literature Number: BMET0210.0 REV101512 & 510(k) Summary K090871.
- 20. Ref. Wright Medical Literature Number: AP-002461B 21-Nov-2018 and 510K Summary K182823.
- 21. Ref. Stryker Literature Number: 90-07900 LOT B1008 & 510(k) Summary K060061.
- 22. Ref. DePuy Synthes Literature Number: J11469-A & 510(k) Summary K102722.
- 23. Ref. http://subchondroplasty.com/healthcare-professionals-indications.html ZimmerBiomet literature number 903.051.13 FDA 510(K) Summary K101557.
- 24. Ref. Skeletal Kinetics Literature Number: LBL 10208 Rev AD & 510(k) Summary K100986.
- 25. Ref. Bone Solutions Test Report: TD-328 [A] OsteoCrete Specification Setting and Sample Size 5cc

Claims based on critically sized rabbit lateral condyle defect model, rabbit anterior cruciate ligament reconstruction, equine metacarpal and metatarsal fracture fixation, and equine metatarsal osteotomy. It is unknown how results from the rabbit or equine models compare with clinical results in humans.



100 Capital Drive Warsaw, IN 46582 +1-574-306-4006

WishBoneMedical.com

Manufactured for: **WishBone Medical, Inc.** 

All trademarks herein are the property of WishBone Medical, Inc. or its subsidiaries unless otherwise indicated. This material is intended for the sole use and benefit of Health Care Professionals and the WishBone Medical Sales Force. It is not to be redistributed, duplicated or disclosed without the express written consent of WishBone Medical. All WishBONE Mg kits are only available as single-use, sterile packed kits. Always confirm product expiration date prior to use.

For product information, including indications, contraindications, warnings, precautions and potential adverse effects, visit WishBone Medical's Instructions for Use page online: www.WishBoneMedical.com/IFU.