

THE FIRST AND ONLY
**CLEFT PALATE
CONVENIENCE KIT**

WishBone Convenience Kits: Sterile Packed, Ready to Go

Your team aims to provide optimal treatment and care to every patient who needs it—when they need it. This is made possible when surgeons and supporting staff are equipped to work efficiently.

Designed to streamline a palatoplasty from start to finish, the WishBone Cleft Palate Convenience Kit comes sterile packed and provides most of the medical and surgical items needed for surgery. By eliminating the 'scavenger hunt' for various components, these convenience kits enable teams to devote time where it matters most.



WISHBONE
MEDICAL

Because *kids* are not just little adults.®

Learn more at WishBoneMedical.com

CLEFT PALATE CORRECTION

A cleft palate occurs if the tissue that makes up the roof of the mouth does not join together completely during pregnancy. For some babies, both the front and back parts of the palate are open. For others, only part of the palate is open.¹

- About 1 in every 1,700 babies is born with cleft palate in the United States.²

Surgery to repair a cleft palate is recommended within the first 18 months of life or earlier if possible.³

Surgical repair can improve the look and appearance of a child's face and might also improve breathing, hearing, and speech and language development.⁴



ORDER INFORMATION

ITEM DESCRIPTION	QTY
Bottle Oral Rinse	1
Iodine Swabs	1
4x4" X-ray Sponge	1
Knife Handle	1
Adson Forcep	1
Needle w/ Syringe 27g x 1 1/2" 3cc	2
Needle Holder	1
Tenotomy Scissors	1
Periosteal Freer	1
Skin Marker	1
Split Pack	1
3 Compartment Tray	1
Skin Hook	2

KIT #: CP-KIT



Manufactured for WishBone Medical, Inc.

1. Centers for Disease Control and Prevention. (2020, December 28). Facts about cleft lip and cleft palate. Centers for Disease Control and Prevention. <https://www.cdc.gov/ncbddd/birthdefects/cleftlip.html>.

2. Mai CT, Isenburg JL, Canfield MA, Meyer RE, Correa A, Alverson CJ, Lupo PJ, Riehle-Colarusso T, Cho SJ, Aggarwal D, Kirby RS. National population-based estimates for major birth defects, 2010–2014. Birth Defects Research. 2019; 111(18): 1420-1435.

3. American Cleft Palate-Craniofacial Association. Parameters for evaluation and treatment of patients with cleft lip/palate or other craniofacial anomalies. Revised edition, Nov 2009. Chapel Hill, NC. P.1-34.

4. D'Antonio, L., Scherer, N. J. (2019, September 20). Language and play development in toddlers with cleft lip and/or palate. ASHA Wire. <https://pubs.asha.org/doi/10.1044/1058-0360.0604.48>.

Convenience kit manufactured by CPT Medical, Inc.
6000A Pelham Rd., Greenville, SC 29615



WISHBONE
MEDICAL

100 Capital Drive, Warsaw, IN 46582

www.WishBoneMedical.com

+1-574-306-4006