


# MALLEOLAR 4.0mm CANNULATED SCREW KITS

Single-Use. Sterile Packed. Ready to Use.

## IMPLANT FEATURES

- Ti-6Al-4V titanium alloy
- Self-tapping flutes
- Thread length = 1/3 screw length
- 36mm–44mm in 2mm increments; 50mm
- Threaded and non-threaded wires available
- Hemispherical head designed to have annular contact with washers
- Low-profile head designed to recess into washer to reduce soft tissue irritation
- Reverse-cutting flutes designed to assist in screw removal



## KIT COMPONENTS

- |  |   |
|--|---|
| 1 4.0 Cannulated Screws, Short Thread (x2) | 5 Parallel Wire Guide                               |
| 2 Washer, 7.0mm (x2)                       | 6 Direct Measurement Device                         |
| 3 1.25mm Guide Wires 150mm Length          | 7 Cannulated T15 Screwdriver                        |
| 4 2.7mm Cannulated Drill Bit               | 8 Cannulated T15 AO Quick Connect Screwdriver Shaft |



Select system components are available as sterile packed ancillary items. For more information, please refer to the WishBone Product Catalog (LIT-PC-WBM) or contact your local rep for product availability.



6mm diameter head with  
T15 hexalobular drive



## ORDER INFORMATION

### Malleolar Procedure Kit, 4.0mm x 36mm Cannulated Screw, Short Thread with Instruments, Ti

ITEM DESCRIPTION	QTY
4.0mm x 36mm Short Thread Cannulated Screw, Ti	2
7.0mm Washer	2
4.0mm Cannulated Screw Parallel Wire Guide	1
4.0mm Cannulated Screw Direct Measurement Device	1
2.7mm Cannulated Drill Bit	1
T15 Cannulated Screwdriver	1
T15 AO Quick Connect Driver Shaft	1
1.25mm x 150mm Smooth Guide Wire	2
1.25mm x 150mm Threaded Guide Wire	2

KIT # MPK-CSST4036T

### Malleolar Procedure Kit, 4.0mm x 38mm Cannulated Screw, Short Thread with Instruments, Ti

ITEM DESCRIPTION	QTY
4.0mm x 38mm Short Thread Cannulated Screw, Ti	2
7.0mm Washer	2
4.0mm Cannulated Screw Parallel Wire Guide	1
4.0mm Cannulated Screw Direct Measurement Device	1
2.7mm Cannulated Drill Bit	1
T15 Cannulated Screwdriver	1
T15 AO Quick Connect Driver Shaft	1
1.25mm x 150mm Smooth Guide Wire	2
1.25mm x 150mm Threaded Guide Wire	2

KIT # MPK-CSST4038T

### Malleolar Procedure Kit, 4.0mm x 40mm Cannulated Screw, Short Thread with Instruments, Ti

ITEM DESCRIPTION	QTY
4.0mm x 40mm Short Thread Cannulated Screw, Ti	2
7.0mm Washer	2
4.0mm Cannulated Screw Parallel Wire Guide	1
4.0mm Cannulated Screw Direct Measurement Device	1
2.7mm Cannulated Drill Bit	1
T15 Cannulated Screwdriver	1
T15 AO Quick Connect Driver Shaft	1
1.25mm x 150mm Smooth Guide Wire	2
1.25mm x 150mm Threaded Guide Wire	2

KIT # MPK-CSST4040T

### Malleolar Procedure Kit, 4.0mm x 42mm Cannulated Screw, Short Thread with Instruments, Ti

ITEM DESCRIPTION	QTY
4.0mm x 42mm Short Thread Cannulated Screw, Ti	2
7.0mm Washer	2
4.0mm Cannulated Screw Parallel Wire Guide	1
4.0mm Cannulated Screw Direct Measurement Device	1
2.7mm Cannulated Drill Bit	1
T15 Cannulated Screwdriver	1
T15 AO Quick Connect Driver Shaft	1
1.25mm x 150mm Smooth Guide Wire	2
1.25mm x 150mm Threaded Guide Wire	2

KIT # MPK-CSST4042T

### Malleolar Procedure Kit, 4.0mm x 44mm Cannulated Screw, Short Thread with Instruments, Ti

ITEM DESCRIPTION	QTY
4.0mm x 44mm Short Thread Cannulated Screw, Ti	2
7.0mm Washer	2
4.0mm Cannulated Screw Parallel Wire Guide	1
4.0mm Cannulated Screw Direct Measurement Device	1
2.7mm Cannulated Drill Bit	1
T15 Cannulated Screwdriver	1
T15 AO Quick Connect Driver Shaft	1
1.25mm x 150mm Smooth Guide Wire	2
1.25mm x 150mm Threaded Guide Wire	2

KIT # MPK-CSST4044T

### Malleolar Procedure Kit, 4.0mm x 50mm Cannulated Screw, Short Thread with Instruments, Ti

ITEM DESCRIPTION	QTY
4.0mm x 50mm Short Thread Cannulated Screw, Ti	2
7.0mm Washer	2
4.0mm Cannulated Screw Parallel Wire Guide	1
4.0mm Cannulated Screw Direct Measurement Device	1
2.7mm Cannulated Drill Bit	1
T15 Cannulated Screwdriver	1
T15 AO Quick Connect Driver Shaft	1
1.25mm x 150mm Smooth Guide Wire	2
1.25mm x 150mm Threaded Guide Wire	2

KIT # MPK-CSST4050T

If "Made in USA" product configuration is required, please add "-US" to kit number when ordering.



100 Capital Drive  
Warsaw, IN 46582  
+1-574-306-4006

For product information, including indications, contraindications, warnings, precautions and potential adverse effects, visit WishBone Medical's Instructions for Use page online: [www.WishBoneMedical.com/IFU](http://www.WishBoneMedical.com/IFU).