

2.7MM & 3.5MM CORTICAL SCREW KITS


Single-Use. Sterile Packed. Ready to Use.

WishBone Medical 2.7mm and 3.5mm Cortical Screws are used for pediatric and adult patients as indicated for pelvic, small and long bone fracture fixation and fixation of bones that have been surgically prepared (osteotomy) for correction of deformity or arthrodesis.

Complete procedure kits in single-use sterile packaging eliminate set processing and help avoid delays or cancellations attributed to missing components.

Readily available, pristine instruments—brand new every time!

- | | |
|----------------------------|-------------------------|
| 1 Cortical Screws* | 4 Small Depth Gauge |
| 2 Double-Ended Drill Guide | 5 Solid T15 Screwdriver |
| 3 2.0mm (or 2.8) Drill Bit | |



*Non-locking screws included in kit. Ancillary screws in 8-50mm lengths available separately.

NON-LOCKING CORTICAL SCREWS

- Designed for fixation in cortical bone
- All screws utilize a T15 driver



Because kids are not just little adults.®

Learn more at WishBoneMedical.com

ORDER INFORMATION

2.7mm Cortical Screws (10mm-16mm) Procedure Kit

ITEM DESCRIPTION	QTY
2.7mm x 10mm Ti Cortical Screw	1
2.7mm x 12mm Ti Cortical Screw	1
2.7mm x 14mm Ti Cortical Screw	1
2.7mm x 16mm Ti Cortical Screw	1
Small Depth Gauge (65mm)	1
2.0mm x 150mm solid drill bit	1
Solid T15 Screwdriver	1
2.0/2.7mm Double-Ended Drill Guide	1

KIT # PK-TCX27SM

2.7mm Cortical Screws (18mm-24mm) Procedure Kit

ITEM DESCRIPTION	QTY
2.7mm x 18mm Ti Cortical Screw	1
2.7mm x 20mm Ti Cortical Screw	1
2.7mm x 22mm Ti Cortical Screw	1
2.7mm x 24mm Ti Cortical Screw	1
Small Depth Gauge (65mm)	1
2.0mm x 150mm solid drill bit	1
Solid T15 Screwdriver	1
2.0/2.7mm Double-Ended Drill Guide	1

KIT # PK-TCX27LG

3.5mm Cortical Screws (10mm-16mm) Procedure Kit

ITEM DESCRIPTION	QTY
3.5mm x 10mm Ti Cortical Screw	1
3.5mm x 12mm Ti Cortical Screw	1
3.5mm x 14mm Ti Cortical Screw	1
3.5mm x 16mm Ti Cortical Screw	1
Small Depth Gauge (65mm)	1
2.8mm x 150mm solid drill bit	1
Solid T15 Screwdriver	1
2.8/3.5mm Double-Ended Drill Guide	1

KIT # PK-TCX35SM

3.5mm Cortical Screws (18mm-24mm) Procedure Kit

ITEM DESCRIPTION	QTY
3.5mm x 18mm Ti Cortical Screw	1
3.5mm x 20mm Ti Cortical Screw	1
3.5mm x 22mm Ti Cortical Screw	1
3.5mm x 24mm Ti Cortical Screw	1
Small Depth Gauge (65mm)	1
2.8mm x 150mm solid drill bit	1
Solid T15 Screwdriver	1
2.8/3.5mm Double-Ended Drill Guide	1

KIT # PK-TCX35LG

Select system components are available as sterile packed ancillary items. For more information, please refer to the WishBone Product Catalog (LIT-PC-WBM) or contact your local rep for product availability.

If "Made in USA" product configuration is required, please add "-US" to part number when ordering.



WishBoneMedical.com

100 Capital Drive
Warsaw, IN 46582
+1-574-306-4006

For product information, including indications, contraindications, warnings, precautions and potential adverse effects, visit WishBone Medical's Instructions for Use page online: www.WishBoneMedical.com/IFU.