

# PEDIATRIC COMBO K-WIRE SYSTEM

## The Single-Use Kit Every OR Needs for K-Wire Pinning Procedures

The WishBone K-Wire System, is intended to be used for fixation of pediatric bone fractures, bone reconstruction; and guide pins for insertion of other implants.

- Two kits for various pinning procedures
- Smooth 316 stainless steel wires are simply cut-to-length @ 200mm, minimizing inventory
- Precision single-ended trocar tips



### SYSTEM COMPONENTS

- 1 Double K-wire Soft Tissue Protector, 1.25mm/1.60mm (or 1.60/2.00mm)
- 2 1.25mm (or 1.60) K-wire, 200mm (x3)
- 3 1.60mm (or 2.00) K-wire, 200mm (x3)
- 4 Hex Driver, 1.5mm
- 5 K-wire Sports Balls, 3 pack (x2)



WISHBONE  
MEDICAL

## EVERYTHING YOU NEED. BRAND NEW EVERY TIME.

WishBone's single-use K-Wire Procedure Kits come sterile packed and ready to use, eliminating the need to sterilize implants & instruments and decreasing overall costs.

Plus, labels are included in each kit for convenient billing and inventory management.



## ORDER INFORMATION

### Elbow Fracture Treatment Procedure Kit

#### KIT NUMBER

PK-KWELB-US	PK-KWELB	DESCRIPTION	QTY
13-904-002020S-1**	13-904-002020S-0**	2.00mm wire, 200mm	3
13-904-001620S-1	13-904-001620S-0	1.60mm wire, 200mm	3
51-904-101620G-1	—	Double Soft Tissue Protector, 1.60mm/2.00mm	1
HDRASS15-US	HDRASS15	Hex Driver for K-Wire Sports Balls Tip Protectors, 1.5mm	1
KWSB-US	KWSB	Sports Balls, 3 Pack	2

US MADE

OUS MADE

\*\*Each package contains 3 wires.

ANCILLARY\*  
ITEM NUMBER

### Extremity Fracture Treatment Procedure Kit

#### KIT NUMBER

PK-KWEXT-US	PK-KWEXT	DESCRIPTION	QTY
13-904-001220S-**	13-904-001220S-0**	1.25mm wire, 200mm	3
13-904-001620S-1	13-904-001620S-0	1.60mm wire, 200mm	3
51-904-101620G-1	—	Double Soft Tissue Protector, 1.25mm/1.60mm	1
HDRASS15-US	HDRASS15	Hex Driver for K-Wire Sports Balls Tip Protectors, 1.5mm	1
KWSB-US	KWSB	Sports Balls, 3 pack	2

US MADE

OUS MADE

\*\*Each package contains 3 wires.

ANCILLARY\*  
ITEM NUMBER

**New!**



#### SPORTS THEMED WIRE BALLS

- All of the benefits of traditional K-wire balls, with a look that kids will love
- **One size fits all** – Hex driver designed to adjust each ball to fit any K-wire within the system

\*Select system components are available as sterile packed ancillary items. For more information, please refer to the WishBone Product Catalog (LIT-PC-WBM) or contact your local rep for product availability.



WishBoneMedical.com

100 Capital Drive  
Warsaw, IN 46582

+1-574-306-4006

For product information, including indications, contraindications, warnings, precautions and potential adverse effects, visit WishBone Medical's Instructions for Use page online: [www.WishBoneMedical.com/IFU](http://www.WishBoneMedical.com/IFU).