

# MALLEOLAR 4.0

CANNULATED SCREW KITS

Low Profile

Reverse-cutting flutes


Self-tapping

Spherical underside

1/3 Partial Thread



T15 Hexalobe



# MALLEOLAR 4.0

CANNULATED SCREW KITS

## PROCEDURE KITS

- 1 4.0 Cannulated Screws 2 of Each  
(36mm, 38mm, 40mm, 42mm, 44mm, 50mm)
- 2 7.0mm Washer 2 of Each
- 3 1.25mm x 150mm Smooth Guide Wire\* 2 of Each
- 4 1.25mm x 150mm Threaded Guide Wire\* 2 of Each
- 5 2.7mm Cannulated Drill Bit\*
- 6 Parallel Wire Guide
- 7 Direct Measurement Device
- 8 Cannulated T15 Screwdriver\*
- 9 Cannulated T15 Driver Shaft, AO\*

\*Available as individually packaged ancillary items



SINGLE-USE.  
STERILE PACKED.  
READY TO USE.



WISHBONE  
MEDICAL

100 Capital Drive  
Warsaw, IN 46582  
+1-574-306-4006  
WishBoneMedical.com



© 2022 WishBone Medical, Inc. | LIT-SB-40CSS REV C | 10-2022

We **WORK** so they can **PLAY**.™

Select system components are available as sterile packed ancillary items.

For ordering information, please contact your local rep or visit the WishBone Medical product catalog page online at [www.wishbonemedical.com/catalog](http://www.wishbonemedical.com/catalog)



For product information, including indications, contraindications, warnings, precautions and potential adverse effects, visit WishBone Medical's Instructions for Use page online: [www.WishBoneMedical.com/IFU](http://www.WishBoneMedical.com/IFU)