

CORTiCAL Screw Kits



All implants made from
medical-grade titanium alloy

CORTICAL Screw Kits

PROCEDURE KITS

2.7mm Kits

- 1 **2.7mm Cortical Screws***
(10, 12, 14, 16mm or 18, 20, 22, 24mm)
- 2 2.0/2.7mm Double-Ended Drill Guide
- 3 2.0mm Drill Bit*
- 4 Small Depth Gauge*
- 5 Solid T15 Screwdriver*

3.5mm Kits

- 1 **3.5mm Cortical Screws***
(10, 12, 14, 16mm or 18, 20, 22, 24mm)
- 2 2.8/3.5mm Double Ended Drill Guide
- 3 2.8mm Drill Bit*
- 4 Small Depth Gauge*
- 5 Solid T15 Screwdriver*

*Available as individually packaged ancillary items

ANCILLARY SCREWS**

2.7mm or 3.5mm Cortical

8-40mm - 2mm Increments

40-50mm - 5mm Increments

**Individually packaged as separate ancillary items



SINGLE-USE.
STERILE PACKED.
READY TO USE.



WISHBONE
MEDICAL

100 Capital Drive
Warsaw, IN 46582
+1-574-306-4006
WishBoneMedical.com



© 2022 WishBone Medical, Inc. | LIT-SB-CSK REV B | 10-2022

We **WORK** so they can **PLAY**™

Select system components are available as sterile packed ancillary items.

For ordering information, please contact your local rep or visit the WishBone Medical product catalog page online at www.wishbonemedical.com/catalog

For product information, including indications, contraindications, warnings, precautions and potential adverse effects, visit WishBone Medical's Instructions for Use page online: www.WishBoneMedical.com/IFU

