

K-Wire

SYSTEM

PEDIATRIC K-WIRE SYSTEM

The Single-Use Kit Every O.R. Needs for K-Wire Pinning Procedures

The WishBone K-Wire System is intended to be used for fixation of pediatric bone fractures, bone reconstruction; and guide pins for insertion of other implants.



FEATURES

- Three K-Wire kits size offering for various pinning procedures
- All wires are made from smooth 316 stainless steel
- Wires are 200mm long to minimize inventory
- Precision single-ended trocar tip
- Eliminates the need to sterilize implants & instruments
- Labels included in each system for convenient billing & inventory management



EVERYTHING YOU NEED. BRAND NEW EVERY TIME.

WishBone's single-use K-Wire Procedure Kits come sterile packed and ready to use, eliminating the need to sterilize implants & instruments and decreasing overall costs.

Plus, labels are included in each kit for convenient billing and inventory management.

Elbow Fracture Kit - 1.60mm

ITEM DESCRIPTION	QTY
1.60mm Wire, 200mm	3
Double Soft Tissue Protector, 1.60mm/2.00mm	1
Hex Driver for K-Wire Sports Balls Tip Protectors, 1.5mm	1
Sports Balls, 3 Pack	1

Elbow Fracture Kit - 2.00mm

ITEM DESCRIPTION	QTY
2.00mm Wire, 200mm	3
Double Soft Tissue Protector, 1.60mm/2.00mm	1
Hex Driver for K-Wire Sports Balls Tip Protectors, 1.5mm	1
Sports Balls, 3 Pack	1

Extremity Fracture Kit - 1.25mm

ITEM DESCRIPTION	QTY
1.25mm K-wire, 200mm	3
Double Soft Tissue Protector, 1.25mm/1.60mm	1
Hex Driver for K-Wire Sports Balls Tip Protectors, 1.5mm	1
Sports Balls, 3 Pack	1

NEW!



SPORTS THEMED WIRE BALLS

- All of the benefits of traditional K-wire balls, with a look that kids will love

SINGLE-USE. STERILE PACKED. READY TO USE.

We **WORK** so they can **PLAY**.™

Select system components are available as sterile packed ancillary items.

For ordering information, please contact your local rep or visit the WishBone Medical product catalog page online at www.wishbonemedical.com/catalog

For product information, including indications, contraindications, warnings, precautions and potential adverse effects, visit WishBone Medical's Instructions for Use page online: www.WishBoneMedical.com/IFU



WISHBONE
MEDICAL

100 Capital Drive
Warsaw, IN 46582
+1-574-306-4006
WishBoneMedical.com

